

**LABORATÓRIOS - V.N.FAMALICÃO**

COMPANY  
**HELIOTEXTIL - ETIQUETAS E PASSAMANARIAS, S.A.**  
**APARTADO 467**

**3702- S. JOÃO DA MADEIRA CODEX**

**Request Nr.: 4722/2011**

**Date of Reception : 2011/07/13**

**Observations**

**Groups of tests**

TOXICIDADE EN 71-PARTE 3

**Samples reference - Your reference**

8420/2011 - "FEROFIX" - Ref: 0823127

**Tests required**

Chemical Tests


**COMMENTS**

The sample is ok.

- The tests were performed between the following dates: 2011/07/13 and 2011/07/18.

V.N.FAMALICÃO, 18th July 2011

LABORATORY  
COORDINATOR



(Deolinda Domingues)

**NOTES:**

- These results were obtained according to the proceedings referred in the Quality Manual of CITEVE and concern only the samples submitted to testing (above mentioned).
- Any part of this report cannot be reproduced without the prior permission of CITEVE.
- The tests signalled by \* are not included in the scope of accreditation of this laboratory
- q.l - quantification limit                      d.l. - detection limit                      n.d. - not detected

LABORATÓRIOS - V.N.FAMALICÃO

<u>Sample Ref.</u>	<u>Your ref.</u>	<u>Sample description</u>
8420/2011	"FEROFIX" - Ref: 0823127	Label sample

**Test/Standard:** \* ARSENIC (TOYS-MIGRATION) / NP EN 71-3:1997 (ICP-MS)

**Results**

Arsenic (mg/kg)	< 0,25 (q.l.)
Analytical method: ICP-MS	

**Test Conditions**

Preparation according to: 8.4

**Test/Standard:** \* ANTIMONY (TOYS-MIGRATION) / NP EN 71-3:1997 (ICP-MS)

**Results**

Antimony (mg/kg)	< 2,5 (q.l.)
Analytical method: ICP-MS	

**Test Conditions**

Preparation according to: 8.4

**Test/Standard:** \* BARIUM (TOYS-MIGRATION) / NP EN 71-3:1997 (ICP-MS)

**Results**

Barium (mg/kg)	< 2,5 (q.l.)
Analytical method: ICP-MS	

**Test Conditions**

Preparation according to: 8.4

**Test/Standard:** \* CADMIUM (TOYS-MIGRATION) / NP EN 71-3:1997 (ICP-MS)

**Results**

Cadmium (mg/kg)	< 0,25 (q.l.)
Analytical method: ICP-MS	

**Test Conditions**

Preparation according to: 8.4

**Test/Standard:** \* LEAD (TOYS-MIGRATION) / NP EN 71-3:1997 (ICP-MS)

**Results**

Lead (mg/kg)	< 0,25 (q.l.)
--------------	---------------

## LABORATÓRIOS - V.N.FAMALICÃO

Analytical method: ICP-MS

### Test Conditions

Preparation according to: 8.4

**Test/Standard:** \* *MERCURY (TOYS-MIGRATION) / NP EN 71-3:1997 (ICP-MS)*

### Results

Mercury (mg/kg) < 0,05 (q.l.)  
Analytical method: ICP-MS

### Test Conditions

Preparation according to: 8.4

**Test/Standard:** \* *SELENIUM (TOYS-MIGRATION) / NP EN 71-3:1997 (ICP-MS)*

### Results

Selenium (mg/kg) < 2,5 (q.l.)  
Analytical method: ICP-MS

### Test Conditions

Preparation according to: 8.4

**Test/Standard:** \* *CHROMIUM (TOYS-MIGRATION) / NP EN 71-3:1997 (ICP-MS)*

### Results

Chromium (mg/kg) < 2,5 (q.l.)  
Analytical method: ICP-MS

### Test Conditions

Preparation according to: 8.4